



Restrictions within CE marking regarding Third-party Trunk support

The restrictions specified below apply to the use of third-party trunk supports while maintaining the CE marking on the following product series from Mercado Medic AB; REAL 9000 PLUS, REAL 6100 PLUS and REAL 6100 PLUS LiNX.

Assembly must be done with Mercado Medic's track system for trunk support mounted on Mercado Medic's backrests intended for this purpose, with M6 thread for the track system.

805575 Track 300mm complete for accessories on backrest (Pcs)

805598 Track 200mm complete for accessories on backrest (Pcs)

Example of backrests for track system:

S-BR1010XX Backrest ErgoMedic 46x43 for track system

- S- = Threaded hole pattern for track system
- XX = Color code: Currently available in the following:
 - 45 = Plush Real Blue
 - 46 = Plush Grey
 - 48 = Plush Black

If the above criteria are met, the CE marking can be maintained on the product. This means that the product can be restored to its original condition without requiring Mercado Medic's approval.

However, the product still has a maximum expected service life of 10 years, after which it requires an inspection according to chapter 6.3 in the Instructions for use for the above-mentioned products.