



THE REAL-CHAIRS ARE CE-CERTIFIED
IN ACCORDANCE WITH THE MPA STATUTES
OF MEDICAL DEVICES, LVFS 2003:11.

PRODUCT RANGE

a selection of the most common REAL-chairs



The Company

Mercado Medic AB is engaged in the business of development, sales and service of special chairs for all people with functional disabilities.

We also produce specially adapted chairs to various clinics, offices and industries.

Our head quarters are situated in Sweden and our export office, which distributes our products to the rest of Europe, in the Netherlands.



The REAL System

The REAL system is a Swedish quality chair with a modular construction for simple individual adjustment.

There is a large variety of seats, backrests, wheels, height range and accessories you can easily combine to create a functional and comfortable chair.

It is suited perfectly to be used at school or home as well as at the workplace.

In our webguide you can build your own chair, find information about different components or accessories and much more .

<http://shop.mercado.se/Welcome/>

The REAL Modular System



REAL 9000 PLUS

Sit in a REAL 9000 PLUS at home and at work. It is easy to adjust position for whatever you want to do.

This model has a gasregulated height adjustment of 20 cm. When the brakes are on, the chair stands firmly in position when you want to get up or transfer. The easily rolling wheels make mobility easy.



REAL 9000 PLUS EL

This model is equipped with electric height adjustment of an impressive 28 cm. By pressing the button you can raise or lower the automatic seat to the desired height.

Brakes and the adjustable angle of the seat and backrest are also available with an electric regulation device.

REAL 9000 PLUS COXIT / EL

This seat is divided with individually adjustable flaps. They are attached without any dust trapping recesses.

This model is best suited for people with various types of orthopaedic injuries.

For chairs with an adjustable seat angle the paddle levers are situated on at the back of the seat. On the Coxit chair the paddle lever is situated on the flap



REAL 9000 PLUS CHILD / EL

Children with various disabilities are able to sit properly at a table with the REAL 9000 PLUS CHILD or REAL 9000 PLUS CHILD EL. The adjustment functions has been developed in consultation with medical experts and provide active functional sitting posture.

The pediatric REAL-chair is available with manual or electric lift as well as a number of functional options and accessories



REAL SPIDER

One of the key features of the REAL spider is its low seat height. This, and the maximum legroom in combination with the smooth-running wheels, allow the child to manoeuvre themselves easily both at home and at school.

The smooth-running wheels allow the child to get to a low desk or table on their own, without their teacher/parent/carer having to push them.

The brake can be operated by hand or by foot. This means that the chair will remain safely in position where you want it, but that it can also easily be moved if required.





REAL 2011

The REAL 2011 is part of a series of standing support chairs for support in a standing working position. By dividing the load between legs, seat and backrest you can avoid injuries caused by high static load. We also provide the REAL 2015 with wheels and brake (as shown in the picture to the right). These chairs can be ordered with a variety of armrests, saddle seats, backrest and seat height.

REAL 2624 / 2625

Fatigue and muscle tension is often caused by sitting leaned forward for longer period of times without the right support. REAL 2624 and REAL 2625 are useful in situations and environments where this position is required. REAL 2624 is built on a 4-legged frame with brake. The REAL 2625 is built on a 5-legged frame and has the option of rotation.

REAL 2624 works as well as a aid for both staff and patient. As a feeding chair both REAL 2624 and REAL 2625 are very functional and flexible. These models are available with a variety of options regarding seat height, saddle seats, stomach support and upholstery to allow for adaptation for various users. The high quality allows for years of use without the need of service or repairs.



REAL 9500/9600 PLUS EL

The REAL 9500 and REAL 9600 are classified as working chairs, as opposed to standing support chairs. They are specifically designed for weaker people that needs to ease the burden on the hips, knees or ankles during chores that require a standing position at home or at work.

These models are built on the smaller PLUS-frame. This makes moving around and working in a smaller space, like a kitchen, alot easier.

The option for a bigger H-fram is available to allow for more room for the feet if the user requires it. Beyond the selection of seats and backrests, there is also the option of a electric brake function.





REAL 6100 PLUS



REAL 9200 TWIN EL

is a module bariatric chair, intended as a manual mobility chair for somewhat heavier users. It is equipped with either double electric- or double gaspring actuator, a double backmechanism and a reinforced X-base.

The electric model, REAL 9200 TWIN EL, fits excellent as a “standing up aid” and has a capacity of up to 275 kg user weight. The manual model, REAL 9200 TWIN, has a capacity of up to 180 kg user weight.

REAL 9000 PLUS EL Extra Low

Since the demand for a good seating solution for short people has increased we now offer REAL 9000 PLUS EL with extra low electric lift. This allows for the seat height of the REAL 9000 PLUS EL to be as low as 32 cm.



REAL 6100 PLUS

is an electrically powered wheelchair designed for indoor use. It is constructed to provide the user maximum freedom of movement at home, at the workplace or at school. The range of electric height adjustment is an impressive 28 cm. This means that one has an extensive freedom of movement even in elevated positions. It has a very narrow chassi with a width of just 57 cm.

REAL 6100 PLUS is available for both children and adults with a number of useful options as well as a great variety of functional accessories.



REAL 9000 PLUS REHAB

The REAL 9000 PLUS REHAB chair is a secure, manageable and easy-to-use medical chair. It eases the rehabilitation process for the patient at a occupational therapy or rehab clinic, whilst providing the staff with a good working environment. The chairs ErgoMedic seat system with the extra tall backrest provides a secure and ergonomic position in the chair.

The chairs push bar makes it easier for the staff to transport the patient in the chair. REAL 9000 PLUS REHAB is equipped with both hand- and footbrake. The armrests are adjustable in both width and height which makes transportation of the patient easier. REAL 9000 PLUS REHAB has hygienic upholstery which is easy to maintain and keep clean. The chair has a capacity of 140 kg. *(The chair in the picture is equipped with a 1/4 tip-up foot ring.)*



REAL 9000 PLUS EL EYE CLINIC

REAL 9000 PLUS EYE CLINIC is a modern electric examination chair, designed with a focus on patient safety, ergonomics, functionality and flexibility. It has a full 28 cm electric lift range governed from the armrest or by using a remote control.

The chair can also be equipped with more electrical functions – electric seat angle, back angle and electric brake.

REAL 8100 X-RAY EL

For patient safety, the real 8100 x-ray EL chair has a convenient, newly designed, fully welded under-frame with a safe, easily operated brake. The brake can be applied using either a pedal or lever, either manually or by remote control. The easy-to-manoeuvre castors, with dustproof SKF bearings in the castors and clamps, make it easier to manoeuvre the chair.

We also have a manual model, real 8000 x-ray, equipped with gaspring height adjustment. *(The chair in the picture is extra equipped with a transparent backrest.)*





WELCOME to our webshop!

Go to www.mercado.se, choose language "English". Then, just click on "WEBSHOP" in the menu.

The screenshot displays the Mercado Medic webshop interface. At the top, there's a navigation bar with the Mercado Medic logo and a language selection menu (Svenska, English, Nederlands, French). A red dotted line highlights the 'WEBSHOP' button in the menu. Below the main website view, a detailed view of the webshop page is shown. The page features a login form on the left with fields for 'Username:' and 'Password:', a 'Login' button, and links for 'Forgotten your username and password?' and 'NEW CUSTOMER?'. A central banner reads 'Welcome to Mercado Medic's shop!'. To the right, there are three main sections: 'Spare parts & Accessories' with a search bar and a chair icon, 'Exploded views' with a search bar and a diagram icon, and 'Easy shopping' with a search bar and a PDF icon. Below these sections, there's a row of images showing chair components: a seat cushion, a chair base, a control panel, and a chair frame. The footer contains contact information for Mercado Medic AB, customer support links (FAQ, Contact Us, Technical Support), a 'REGISTERED' badge (Register, Confidentiality), and 'ADDITIONALLY' links (Easy Shopping, Terms & Conditions, About Mercado Shop). Payment logos for VISA, MasterCard, and DIBS are also present. At the bottom, it states 'YOUR ORDER CONTAINS 0 items' and provides the company address and contact details.



Mercado Medic AB
Tryffelslingan 14
181 57 Lidingö
Sweden

Telephone & Fax
Tel 08 555 143 00
Fax 08 555 143 99

Email & Internet
Email info@mercado.se
Internet www.mercado.se