

Standing support chairs with backrest



REAL 2001



REAL 2011



REAL 2014



REAL 2015

DATA STANDING SUPPORT CHAIRS

*standard

Model	2001	2001-1	2011	2011-1	2014	2014-1	2015	2015-1
Lyftomat.....	Yes*	Yes*	Yes*	Yes*	Yes*	Yes*	Yes*	Yes*
Manual seat angle adjustment (+22°/-14°).....	-	-	Yes*	Yes*	Yes*	Yes*	Yes*	Yes*
Rotation lock.....	-	-	Yes*	Yes*	-	-	-	-
Backrest 33x16 cm.....	Yes*	Yes*	-	-	-	-	-	-
Backrest 33x26 cm.....	-	-	Yes*	Yes*	Yes*	Yes*	Yes*	Yes*
lide buttons.....	Yes*	Yes*	Yes*	Yes*	-	-	-	-
Castor shaped gliding feet 50 mm.....	Yes	Yes	Yes	Yes	-	-	-	-
Weight breakable castors 50 mm.....	Yes	Yes	Yes	Yes	-	-	-	-
Castor 75 mm.....	-	-	-	-	Yes*	Yes*	Yes*	Yes*
Castor 50 mm with valve lock.....	Yes	Yes	Yes	Yes	-	-	-	-
Brake.....	-	-	-	-	Footbreak*	Footbreak*	Handbreak*	Handbreak*
Seat height cm.....	54-74	64-84	60-80	70-90	60-80	70-90	60-80	70-90
Seat round 35 cm.....	Yes*	Yes*	-	-	-	-	-	-
Seat round 40 cm.....	Yes	Yes	-	-	-	-	-	-
Sadle seat 36x36 cm.....	-	-	Yes*	Yes*	Yes*	Yes*	Yes*	Yes*
Alternative sadle seat 41x41 cm.....	-	-	Yes	Yes	Yes	Yes	Yes	Yes
Alternative sadle seat Medic 40x36 cm.....	-	-	Yes	Yes	Yes	Yes	Yes	Yes
Alternative sadle seat Ergomedic 40x36 cm.....	-	-	Yes	Yes	Yes	Yes	Yes	Yes
Base cm (w x d).....	42,5x42,5	42,5x42,5	42,5x42,5	42,5x42,5	43x43	43x43	43x43	43x43
Max user weight.....	130 kg	130 kg	130 kg	130 kg	130 kg	130 kg	130 kg	130 kg
Upholstery								
Plush: grey*, red, bown, green, marine, blue, turquoise								
Skai: red, grey, blue, brown, black								
Cotton manchester: red, grey								
Cover for seat and backrest.....	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes

Stomach support chairs



REAL 2621

Stomach support, 27x23 cm



REAL 2624

Stomach support, 27x23 cm



REAL 2625

Sadle seat 36x36cm with alternative
Stomach support 33x16 cm

DATA STOMACH SUPPORT CHAIRS

*standard

Model	2621	2621-1	2624	2624-1	2625
Lyftomat.....	Yes*	Yes*	Yes*	Yes*	Yes*
Fixed lyftomat.....	-	-	Yes	Yes	-
Manual seat angle adjustment (+22°/-14°).....	Yes*	Yes*	Yes*	Yes*	Yes*
Rotationlock.....	Yes*	Yes*	-	-	Yes
Stomach support Medic 27x23 cm.....	Yes*	Yes*	Yes*	Yes*	Yes
Stomach support low 33x16 cm.....	Yes	Yes	Yes	Yes	Yes
Glide buttons.....	Yes*	Yes*	-	-	-
Castor 75 mm.....	-	-	Yes*	Yes*	Yes
Castor shaped feet 50 mm.....	Yes	Yes	-	-	-
Weight breakable castors 50 mm.....	Yes	Yes	-	-	-
Castor 50 mm with valve lock.....	Yes	Yes	-	-	-
Brake.....	-	-	Fotbroms*	Fotbroms*	-
Seat height cm.....	60-80	70-90	60-80	70-90	50-70
Sadle seat Medic 40x36 cm.....	Yes*	Yes*	Yes*	Yes*	Yes*
Alternative sadle seat Ergomedic 40x36 cm.....	Yes	Yes	Yes	Yes	Yes
Alternative sadle seat 36x36 cm.....	Yes	Yes	Yes	Yes	Yes
Alternative sadle seat 41x41 cm.....	Yes	Yes	Yes	Yes	Yes
Base cm (w x d).....	42,5x42,5	42,5x42,5	43x48	43x48	-
Max user weight.....	130 kg	130 kg	130 kg	130 kg	130 kg

Upholstery

Plush: grey*, red, brown, green, marine, blue, turquoise

Skai: red, grey, blue, brown, black

Cotton manchester: red, grey

Cover for seat and backrest..... Yes

Yes

Yes

Yes

Yes