



REAL EYE PATIENT CHAIRS

Swedish-made electric examination chairs



Valid from 090401

Last change 130729

Article number

REAL 9100 EL CLINIC



- Electric height-/lowering w remote control (more alt. avail.)
- Seat height 47-75 cm (alt. 40-60 cm)
- X-base "clinic" 48x53 cm with built-in foot brake (alt. avail.)
- Wheel 100 mm XL (alt. avail.)
- Without seat tilt (alt. seat tilt avail.)
- Back rest 38x26 cm (alt. avail.)
- Seat 44x48 cm (alt. avail.)
- Armrest 35 cm (alt. avail.)
- Clothing in skai (alt. avail.)



base 9000 "Clinic" 48x53 cm

STABLE FOR THE PATIENT AND EASY TO USE FOR STAFF

In a patient chair high standards are set. It must meet the needs of doctors, nurses and patients.

REAL 9100 EL EYE CLINIC CHAIR is a Swedish-made electric examination chair, designed in view of ergonomics and function.

In a patient chair high standards are set. It must meet the needs of doctors, nurses and patients. The design should be stable and the wheels to roll easily. To simplify the examination the chair should also be easily height adjustable. High quality guarantees long life with minimal need for maintenance is also important.

OPERATE THE CHAIR FROM A DISTANCE

real 9100 el eye clinic chair is a Swedish-made electric examination chair, designed in view of ergonomics, function and flexibility. It is comfortable, stable and easily maintained. The chair has a full 28 cm electric lift range governed from the armrest or by using a remote control.

A remote control allows the examiner to operate the chair from a distance.

It can also be equipped with more electrical functions – electric seat angle, back angle and electric brake.

The chair has a welded X-base "Clinic" with built in foot brake (see picture above).

OPTIONS

REAL 9100 EL EYE CLINIC CHAIR is available with a wide range of accessories. One can choose between different types of seat and back angles, alternative seat and backrest sizes and a large variety of accessories.

REAL 9100 PLUS EL CLINIC

- Electric height-/lowering w remote control (more alt. avail.)
- Seat height 47-75 cm (alt. 40-60 cm)
- PLUS-base 48x53 cm with built-in foot brake (alt. avail.)
- Wheel 100 mm XL (alt. avail.)
- Without seat tilt (alt. seat tilt avail.)
- Back rest 38x26 cm (chair on image has backrest 40x34 cm)
- Seat 44x48 cm (alt. avail.)
- Armrest 35 cm (alt. avail.)
- Clothing in skai (alt. avail.)



narrower PLUS-base 43x53 cm



TWO PLUS-BASES INCREASES THE FLEXIBILITY

Mercado Medic AB has advanced our renowned eye clinic chair. To increase the flexibility, we designed an examination chair where you as a customer can choose between two modern Plus-bases in different sizes depending on need.

REAL 9100 PLUS EL EYE CLINIC CHAIR is a modern Swedish-made electric examination chair, designed with a focus on patient safety, ergonomics, functionality and flexibility.

We put high demands in our new eye clinic chair. It must be able to meet the different needs of each individual eye clinic. The design should be secure but also convenient and flexible. To simplify the examination, the chair should also be easily height adjustable. Long life with minimal need for maintenance is also important.

NEW STABLE PLUS-BASE

REAL 9100 PLUS EL EYE CLINIC CHAIR is a Swedish-made electric examination chair, designed in view of ergonomics, function and flexibility. We have focused on your work environment. Therefore, we have made the chair safe for the patient, but also convenient and easy to care for you. REAL 9100 PLUS EL EYE CLINIC CHAIR has a full 28 cm electric lift range governed from the armrest or by using a remote control. A remote control allows the examiner to operate the chair

from a distance. It can also be equipped with more electrical functions – electric seat angle, back angle and electric brake. The chair has a welded PLUS-base with built in foot brake. It offers the patient more room for his/her feet. The brake pedal means that you, the user, can easily brake the chair without using the brake lever.

A NARROWER PLUS-BASE

In the eye clinics that have narrow examination tables a 5 cm narrower plus base is an option (see image above).

OPTIONS

REAL 9100 PLUS EL EYE CLINIC CHAIR is available with a wide range of accessories. One can choose between different types of seat and back angles, alternative seat and backrest sizes and a large variety of accessories that makes your job a little easier.

REAL EYE PATIENT CHAIRS

a selection of accessories



SELECTION OF ACCESSORIES

- 1 Remote controls, multi (2-4 funktion)
up/down (1-funktion)
- 2 Electric control in separate box (4-funktion)
- 3 Foldable armrest
- 4 Neck support cushion ErgoMedic Plus
- 5 Standard wheels, XL 100 mm
- 6 ESD wheels, 100 mm Antistatic
- 7 Leg support depth 100 mm
(follows height adjustment)
- 8 Footring 1/4, foldable, X-base
Footring 1/4, foldable, Plus-base 48
Footring 1/4, foldable, Plus-base 43
- 9 Armrest ErgoMedic 25, 30, 35
- 10 Base 9000 Plus 43 100 XL
- 11 Base 9000 Plus Clinic 48, 100 XL, White
- 12 Base 9000 X Clinic w built in foot brake, 100 XL
- 13 Extra high backrest ErgoMedic 42x43 cm
- 14 Push bar clinic



Mercado Medic AB
Tryffelslingan 14
181 57 Lidingö
Sweden

Telefon & Fax
Tel 08 555 143 00
Fax 08 555 143 99

Email & Internet
Email info@mercado.se
Internet www.mercado.se