



THE REAL-CHAIRS ARE CE-CERTIFIED
IN ACCORDANCE WITH THE MPA STATUTES
OF MEDICAL DEVICES, LVFS 2003:11.

REAL EYE CLINIC CHAIRS

*Electrically controlled examination chairs
made in our own factories in Sweden*



Practical, ergonomic and efficient eye clinic chairs, both for staff and patients.

REAL 9100 PLUS EL eye clinic chair must meet the majority of needs that an eye clinic has. The chair is flexible and functional, while at the same time being stable and safe. REAL 9100 PLUS EL also offers a long lifetime with a minimum need for maintenance.

Ergonomics, safety and function

We have focused on your working environment. That is why we have made the chair safe for patients, while at the same time flexible and easy to use for you. REAL 9100 PLUS EL EASY and REAL 9100 PLUS EL EFFICIENT have a 28 cm electrical lifting range that is controlled from the arm rest or by a remote control. The remote control gives the physician/nurse the opportunity to control the chair, even at distance. REAL 9100 PLUS EFFICIENT also offers electric seat tilt (it can be tilted both backwards and forwards). The chair can be fitted with more electrical functions, (e.g. electric brake) in order to make examinations more efficient. Everything can be controlled from the remote control.

Avoid conflicts with the market's narrowest clinic base

Our narrowest PLUS base is only 43 cm wide and 53 cm deep. It is the only clinic base on the market that is narrow enough to manage all types of examination tables, whatever their width. Patient chairs with a bulky base can result in conflict when the chair's front base leg is too wide to come in between the base legs on the examination table and close to the examination equipment. The foot pedal at the rear makes it easy to brake the chair without having to use the brake lever.

Accessories

REAL 9100 PLUS EL can be complemented with a wide range of accessories to make your work easier (see the rear of this brochure).



REAL 9100 PLUS EL EASY



[The chair in the photo is fitted with PLUS-base clinic 48x53cm]

- Electrical height adjustment by remote control (cordless or with cord)
- Seat height 47-75 cm (or 40-60 cm)
- Without seat tilt (others available)
- Seat ErgoMedic 44x48 cm (others available)
- Back ErgoMedic 38x26 cm (the illustrated chair has backrest 40x34 cm)
- PLUS-base 43x53 cm (or 48x53 cm)
- Wheels 100mm XL (others available)
- Arm rest 35 cm (others available)
- Skai upholstery cover (other colours available)

REAL 9100 PLUS EL EFFICIENT



[The chair in the photo is fitted with PLUS-base clinic 43x53cm]

- Electrical height adjustment by remote control (cordless or with cord)
- Seat height 47-75 cm (or 40-60 cm)
- Electrically controlled seat tilt (others available)
- Seat ErgoMedic 44x48 cm (others available)
- Back ErgoMedic 38x26 cm (the illustrated chair has backrest 42x43 cm)
- PLUS-base 43x53 cm with foot brake (or 48x53 cm)
- Foot and handbrake (or electric brake)
- Wheels 100mm XL (others available)
- Foldable arm rest 35 cm (others available)
- Skai upholstery cover (other colours available)
- ¼ up-foldable foot ring (option)

REAL 9100 PLUS EL eye clinic chair is a modern Swedish-made electrically-controlled patient chair, designed with a focus on patient safety, ergonomics and functionality. In order to increase flexibility, we have created an examination chair in which our customers can choose between two modern PLUS bases with different widths. In order to increase the efficiency of examinations, the REAL 9100 PLUS EL EFFICIENT is fitted with more electrical functions than the REAL 9100 PLUS EL EASY. The chairs can be individually customised with alternative functions for different needs.

REAL EYE CLINIC CHAIRS – a range of accessories



- 1 Remote controls, multi (2-4 functions) up/down (Single function)
- 2 Remote control with cord (4 function)
- 3 Foldable arm rest
- 4 Sampling arm rest
- 5 SL wheels, extra easy rolling wheels, 100 mm
- 6 BXL wheels, extra broad and easy rolling wheels, 100 mm
- 7 Foot rest 100 mm deep (that follows height regulation)
- 8 Foot ring 1/4, foldable, PLUS 48
Foot ring 1/4, foldable, PLUS 43
- 9 Arm rest plates ErgoMedic 25, 30, 35
- 10 Base 9000 PLUS 43 100 XL
- 11 Base 9000 Plus Klinik 48 100 XL White
- 12 Extra high back ErgoMedic 42x43 cm
- 13 Push handle 'Clinic'

Colour options skai



[31 red]



[32 brown]



[35 blue]



[36 grey]



[38 black]

Would you like another colour?
Contact us!

