

Armrest mechanism

REAL 9000 PLUS, REAL 6100 PLUS

Art.nr: 804538
 Changed: 1605
 Valid from: 0706

Article	Art.nr.
ARM1011	Armrest mech. Plus 32 cm R
ARM1012	Armrest mech. Plus 32 cm L
ARM1111	Armrest mech. Plus 38 cm R
ARM1112	Armrest mech. Plus 38 cm L
ARM1211	Armrest mech. Plus mini 21 cm R
ARM1212	Armrest mech. Plus mini 21 cm L
ARM1311	Armrest mech. Plus mini 32 cm R
ARM1312	Armrest mech. Plus mini 32 cm L
ARM1411	Armrest mech. Plus 6100 32 cm R
ARM1412	Armrest mech. Plus 6100 32 cm L
ARM1511	Armrest mech. Plus 6100 38 cm R
ARM1512	Armrest mech. Plus 6100 38 cm L
ARM7111	Armrest mech. Plus 45 cm R
ARM7112	Armrest mech. Plus 45 cm L
ARM3011	Armrest mech. Foldable 32 cm R
ARM3012	Armrest mech. Foldable 32 cm L
ARM3111	Armrest mech. Foldable 38 cm R
ARM3112	Armrest mech. Foldable 38 cm L
Necessary tools	
Allen key 5 mm	800897

Instructions

Assemble as seen below. Adjust the armrest mechanism in width so that the armrest is close to the body, and in height so that the armrest provides support for the elbow with the shoulder in a relaxed position.

