

Article	Art.nr.
Armrest ErgoMedic	AR10xxx
Armrest ErgoMedic Mini	AR11xxx
Armrest ErgoMedic f foldable	AR15xxx
Armrest ErgoMedic Plus	AR20xxx
Armrest Mobile	AR30xxx
Armrest Mobile 9000 Drilled	AR31xxx
Armrest Hemiplegi 45x12cm	AR4060xx
Armrest Comfort	AR50xxx
Armrest extra soft	AR51xxx
K6S M6x25 black chromated furniture screw M6x25	802039

xx: Decides size(25-40cm) and side(ErgoMedic PLUS)

xx: Decides color, color codes can be found in the "Seat systems" folder.

Necessary tools	Art.nr.
Allen key 5 mm	800897

## Instructions

To mount a new armrest, release the screws on the bottom of the armrest(as shown below) using an Allen key. Mount the new armrest and tighten the screws.

