

# Backrest mechanism

REAL 9000 PLUS, REAL 6100 PLUS, REAL 2000

Art.Nr: 803312  
 Changed: 1701  
 Valid from: 0706

## Article

Std backrest mech. low device r back  
 Medic backrest mech low De155/R240 mm Device R back  
*More mechanisms. are available, for complete list, see "Exploded Diagrams"*

## Art.nr.

BRM1000  
 BRM1600

## Necessary tools

Allen key 5 mm

800897

## Instructions

To mount the backrest mechanism, screw it to the seat plate with four M6 screws(1).  
 To mount the backrest, see assembly instructions "Back-/Stomachrest".

