

REAL 9000 PLUS

Try out/order form for CHAIRS
Accessories are available on separate sheet

| | | | | | | | | |
|---|---|---|---|---|---|---|---|------------|
| A | B | C | D | E | F | G | H | Upholstery |
| 9 | | | | | | | | |

Quote request Order

Belongs to order.....



| | | | |
|------------------|--------|-------------------------|------|
| Delivery address | | Invoice address | |
| Tried out by | | Tried out for/markings | |
| Telephone | E-mail | Purchasers order number | Date |

| Basic equipment – Always fill in! | | ALTERNATIVES | |
|--|--|---|--|
| AB Manual Models 90 <input type="checkbox"/> REAL 9000 PLUS 93 <input type="checkbox"/> REAL 9000 PLUS Child** 97 <input type="checkbox"/> REAL 9000 PLUS Coxit <i>*Battery charger included</i> | | AB Electric Models 91 <input type="checkbox"/> REAL 9000 PLUS EL* 94 <input type="checkbox"/> REAL 9000 PLUS Child EL** 98 <input type="checkbox"/> REAL 9000 PLUS Coxit EL* <i>**PLUS-base 43x53 cm</i> | |
| C Seat angle adjustment (+15° Forward, -8° Backward) 0 <input type="checkbox"/> Without adjustment 2 <input type="checkbox"/> Crank adjustment 3 <input type="checkbox"/> Gas spring adjustment at seat | | Upholstery Plush: <input type="checkbox"/> 41 bordeaux <input type="checkbox"/> 42 brown <input type="checkbox"/> 43 green <input type="checkbox"/> 44 marine <input type="checkbox"/> 45 REAL-blue <input type="checkbox"/> 46 grey <input type="checkbox"/> 48 black Patterned fabric Adult: <input type="checkbox"/> 75 blue Patterned fabric Child: <input type="checkbox"/> 81 red <input type="checkbox"/> 85 blue <input type="checkbox"/> 89 yellow Skai: <input type="checkbox"/> 31 red <input type="checkbox"/> 32 brown <input type="checkbox"/> 35 blue <input type="checkbox"/> 36 grey <input type="checkbox"/> 38 black Incontinence upholstery (add cover) <input type="checkbox"/> Seat <input type="checkbox"/> Backrest <input type="checkbox"/> Leg divider Covers: <input type="checkbox"/> Seat <input type="checkbox"/> Backrest <input type="checkbox"/> Others according to below Alternative cotton covers: <input type="checkbox"/> Child/yellow <input type="checkbox"/> Child/blue <input type="checkbox"/> Incontinence cover, grey | |
| D Low back rest 0 <input type="checkbox"/> Standard 1 <input type="checkbox"/> Medic 2 <input type="checkbox"/> Comfort* 3 <input type="checkbox"/> Electric* (EL) | | Alternative Low back rest WidthxHeight <input type="checkbox"/> 35x20 cm Ergomedic PLUS <input type="checkbox"/> 40x20 cm Ergomedic PLUS Alternative High back rest WidthxHeight <input type="checkbox"/> 42x43 cm Ergomedic <input type="checkbox"/> 30x34 cm Ergomedic PLUS <input type="checkbox"/> 35x43 cm Ergomedic PLUS <input type="checkbox"/> 42x43 cm Ergomedic PLUS <input type="checkbox"/> 46x43 cm Ergomedic PLUS <input type="checkbox"/> 30x34 cm Ergomedic PLUS with wedge <input type="checkbox"/> 35x43 cm Ergomedic PLUS with wedge <input type="checkbox"/> 42x43 cm Ergomedic PLUS with wedge <input type="checkbox"/> 46x43 cm Ergomedic PLUS with wedge | |
| D Large back rest 5 <input type="checkbox"/> Standard 6 <input type="checkbox"/> Medic <input type="checkbox"/> Double 7 <input type="checkbox"/> Comfort* 8 <input type="checkbox"/> Electric* (EL) <i>*only for seat 43 cm depth or more</i> | | Alternative seat and back rest <input type="checkbox"/> Seat, only wooden plate <input type="checkbox"/> Seat, only wooden plate with anti glide <input type="checkbox"/> Seat with slow recovery <input type="checkbox"/> Back rest with slow recovery <input type="checkbox"/> All seat angle adj. backwards (0-23°)* <input type="checkbox"/> Adj. seat angle (+15° Front, -23° Back) <input type="checkbox"/> Reversed angle adj. (+8° Front, -15° Back) <i>*Rear base extender is automatically assembled</i> | |
| E Castors 0 <input type="checkbox"/> Standard XL 100 mm 2 <input type="checkbox"/> Extra easy rolling, SL 100 mm* 3 <input type="checkbox"/> Extra wide castors 100 mm <i>*For hard surfaces</i> | | 4 <input type="checkbox"/> 100 mm ESD (anti-static) 5 <input type="checkbox"/> 125 mm** 6 <input type="checkbox"/> 150 mm** <i>**Rear base extender is automatically assembled</i> | |
| F Seat height* Manual EL 0 <input type="checkbox"/> Standard 40-60 cm 43-71 cm 1 <input type="checkbox"/> Low 32-46 cm 36-56 cm 2 <input type="checkbox"/> Extra low 32-47 cm | | 5 <input type="checkbox"/> Standard seat height with rotation 6 <input type="checkbox"/> Low seat height with rotation <i>* Measured without seat angle from floor to underside of seat, with seat in flat position. With seat tilt +1 cm. With seat tilt "adjustable reinforced seat frame" +2 cm.</i> | |
| G Seat ErgoMedic | | Width x Depth Recommended Backrests Armrest | |
| 0 <input type="checkbox"/> Adult <input type="checkbox"/> Coxit 44x48 cm 1 <input type="checkbox"/> Adult <input type="checkbox"/> Coxit 48x48 cm 2 <input type="checkbox"/> Adult <input type="checkbox"/> Coxit 44x43 cm 3 <input type="checkbox"/> Adult/child <input type="checkbox"/> Coxit 40x43 cm 4 <input type="checkbox"/> Adult/child 40x39 cm 5 <input type="checkbox"/> Child (7-12 years) 36x39 cm* 6 <input type="checkbox"/> Child (3-7 years) 33x35 cm* 7 <input type="checkbox"/> Child (7-15 years) 36x43 cm* 8 <input type="checkbox"/> Child 29x32 cm** 9 <input type="checkbox"/> Adult/child <input type="checkbox"/> Coxit 40x48 cm 9 <input type="checkbox"/> Adult 44x53 cm 9 <input type="checkbox"/> Adult 48x53 cm <input type="checkbox"/> Seat in ErgoMedic PLUS design | | Low 38x26 cm high 40x34 cm 35 cm Low 38x26 cm high 40x34 cm 35 cm Low 38x26 cm high 40x34 cm 35 cm Low 33x26 cm high 33x36 cm 35 cm Low 33x26 cm high 33x36 cm 30 cm Low 30x26 cm high 30x34 cm 30 cm Low 27x22 cm high 27x32 cm 25 cm Low 30x26 cm high 30x34 cm 30 cm Low 23x20 cm high 27x32 cm 25 cm Low 33x26 cm high 33x36 cm 35 cm Low 38x26 cm high 40x34 cm 35 cm Low 38x26 cm high 40x34 cm 35 cm <i>* Works with PLUS-base 43x53 cm ** Works with PLUS-base Mini 43x48 cm</i> | |
| H User's weight 0 <input type="checkbox"/> 55-90 kg Manual/EL-lyftomat 1 <input type="checkbox"/> < 55 kg Manual/EL-lyftomat <i>User's weight>140kg contact Mercado</i> | | 2 <input type="checkbox"/> 90-140 kg Manual fixlyftomat 3 <input type="checkbox"/> 90-150 kg EL-lyftomat | |
| EXTRA | | Number of Chairs | |

REAL 9000 PLUS/1409

