



THE REAL-CHAIRS ARE CE-CERTIFIED  
IN ACCORDANCE WITH THE MPA STATUTES  
OF MEDICAL DEVICES LVFS 2003:11.

# REAL 9500/9600 EL

*relieve hips, knees and ankles*



REAL 9500/9600 EL  
is available with PLUS-  
stand 43 cm.

(Note! Not standard.)



## SPECIFICATION

- Seat height adjustment placed under the armrest
- Crank seat angle adjustment
- Standard back rest mechanism (alt. Medic)
- H-frame (alt. PLUS-frame 43 cm)
- XL-casters 100 mm (alt. SL-, BXL-casters 100 mm)
- Seat height 9500: 50-70 cm (alt. low 45-65 cm)
- Seat height 9600 EL: 54-82 cm (alt. low 47-75 cm)
- Saddle seat ErgoMedic 40x36 cm (alt. 36x36 cm)
- Back rest ErgoMedic 33x26 cm (alt. 33x36 cm)
- Back rest ErgoMedic PLUS 35x20 cm (extra option)
- Electric brake on 9600 EL
- Armrest standard: 35 cm (alt. 30, 25 cm)
- Brake lever standard: 45 cm (alt. 35 cm)
- Clothing: plush, imitation leather

## Standing and sitting support

When the muscles are weak, you are in need of a chair that provides good relief. Great demands on function and quality is then necessary.

When your legs are not strong enough, you need a chair both for sitting in and for getting around. This places great demands on function and quality.

The chair has to be comfortable with a sitting position that is ergonomically correct. The wheels should roll easily for easy mobility with reliable brakes which hold the chair firmly in place when you want to stand up. At the same time the chair's functions should be easily adjustable.

Basically, the chair has to work well in every situation.

### *Especially adapted*

REAL 9500/9600 EL are chairs that are especially adapted for people with weak muscles who need to relieve their hips, knees or ankles while working in a standing position at home or at work. These chairs are equipped with, among other things, a H-frame, saddle seat as well as a lifted seat position.



Mercado Medic AB  
Tryffelslingan 14  
181 57 Lidingö  
Sweden

Telefon & Fax  
Tel +46 (0)8 555 143 00  
Fax +46 (0)8 555 143 99

Email & Internet  
Email [info@mercado.se](mailto:info@mercado.se)  
Internet [www.mercado.se](http://www.mercado.se)