



THE REAL-CHAIRS ARE CE-CERTIFIED
IN ACCORDANCE WITH THE MPA STATUTES
OF MEDICAL DEVICES LVFS 2003:11.

REAL[®] 2000

standing support chairs



REAL[®] 2011



- Manual fix liftomat
- Manual angle adjusted seat (prepared for armrest)
- Saddle seat 36x36 cm, Ergomedic
- Backrest 33x26 cm, Ergomedic
- Backrest mechanism Standard
- Seat height 50-70 cm (opt. 60-80 cm)
- User weight 0-130 kg
- Sliding buttons
- Clothing in grey plush

An adaptable standing support chair

For persons that live with a condition of pain or an impaired muscular system, a stand-supporting chair is a good solution. The REAL 2011 is an adaptable stand-supporting chair that makes it possible to have a comfortable half-standing sitting posture.

For persons that live with a condition of pain or an impaired muscular system a stand-supporting chair offers relief with increased function and well-being and the possibility to have a half-standing working posture with the support of a chair. Because the bodyweight is evenly distributed over legs, sitting area and backrest an unnecessary pressure can be avoided.

UNLOAD WITH REAL 2011

The REAL 2011 is an adaptable stand-supporting chair that makes it possible to have a comfortable sitting posture in relation to the work task at hand.

To perform jobs by the kitchen sink or the ironing board are examples on two tasks that require a stand up working posture. Both of which take a long time to fulfill and tiredness in legs and back can become very disturbing. In cases like this a



[Prepared for armrest.]

REAL 2011 can offer the right relief. The chair has a gas-spring regulated interval of 20 cm (70-90 cm alt. 60-80 cm) and a manual angle adjusted seat.

One can easily lock the rotation of the seat with a grip. The REAL 2011 is prepared for mounting of armrests and can be delivered with a wide range of accessories and choices. One can choose between different types of saddle seats, backrests, armrests and a big amount of accessories.

REAL[®] 2015

- Gas-spring controlled liftomat
- Manual angle adjusted seat
- Saddle seat 36x36 cm (opt. available)
- Backrest 33x26 cm (opt. available)
- Seat height 50-70 (opt. 60-80 cm)
- Handbrake
- Frame 43x48 cm
- Castors 75 mm
- Clothing in grey plush



Relief when working in a standing position

REAL 2015 is a stand-supporting chair with castors and handbrake for support when working in a standing position.

This means that you have the possibility to stand with support in your everyday life.

Because the bodyweight is evenly distributed over legs, sitting area and backrest, an unnecessary pressure can be avoided.

Possible to style the workplace

The REAL 2015 stand-supporting chair makes it possible to style the workplace after a comfortable movement pattern.

The shaping of the workplace, regarding reach, disposition and ergonomics can easily be adapted with the half standing sitting posture. The REAL 2015 offers a stable benchmark for the task.

Dentists, photographers, machine operators, servicemen, laboratorie personnel and physiotherapists are examples on users of the REAL 2015.

Support with increased function

To perform jobs by the kitchen sink or the ironing board are examples on two tasks that require a stand up working posture. Both of which take a long time to fulfill and tiredness in legs and back can become very disturbing.

For persons that live with a condition of pain or an impaired muscular system, the REAL 2015 offers relief with increased function and well-being.

REAL[®] 2000



REAL 2001

Max. user weight	130 kg
Seat tilt	–
Rotation	Yes
Seat height	High 60-80 cm
Seat height	Low 50-70 cm
Sliding buttons	Yes
Round seat 35 cm	Yes
Brake	–

REAL 2011

Max. user weight	130 kg
Seat tilt	Ja
Rotation	Ja
Seat height	High 60-80 cm
Seat height	Low 50-70 cm
Sliding buttons	Yes
Saddle seat 36x36 cm	Yes
Brake	–

REAL 2015

Max. user weight	130 kg
Seat tilt	Yes
Rotation	–
Seat height	High 60-80 cm
Seat height	Low 50-70 cm
Castors 75 cm	Yes
Saddle seat 36x36 cm	Yes
Handbrake	Yes

