

REAL[®] 9500/9600 EL

When your legs are not strong enough, you need a chair both for sitting in and for getting around. This places great demands on function and quality. The chair has to be comfortable with a sitting position that is ergonomically correct. The wheels should roll easily for easy mobility with reliable brakes which hold the chair firmly in place when you want to stand up. At the same time the chair's functions should be easily adjustable. Basically, the chair has to work well in every situation.



REAL 9500 & 9600 EL are the chairs that are especially adapted for people with weak muscles who need to relieve their hips, knees or ankles while working in a standing position at home or at work. These chairs are equipped with, among other things, a H-frame, saddle seat as well as a lifted seat position.

REAL 9500/9600 EL

Technical data

- Both REAL 9500 and 9600 EL have the seat height adjustment placed under the armrest as standard
- The chair is delivered with crank seat angle adjustment as standard
- Standard- and medic back rest mechanism is available as alternatives
- H-frame is equipped with either XL- or SL-casters 100 mm
- H-frame with extended rear reduces risk of tipping over
- Seat height:
 - Manual: standard seat height (60-80 cm), alt. low (55-71 cm)
 - Electric: standard seat height (62-90 cm), alt. low (55-75 cm)
- Seat standard: saddle seat ERGOMEDIC 40x36 cm (alt. 36x36 cm, 41x41 cm)
- Back rest 33x26 cm (alt. 33x36 cm)
- ERGOMEDIC
- Arm rest standard: 35 cm (alt. 30, 25 cm)
- Brake lever standard: 45 cm (alt. 35 cm)



Mercado Medic AB
Box 1074
181 22 Lidingö
Sweden

Tel: +46-(0)-8 555 143 00
Fax: +46-(0)-8 555 143 99
E-mail: info@mercado.se

Mercado Medic Export Office:
Nobelweg 22
3899 BN Zeewolde
The Netherlands

Tel: +31-(0)-36 521 99 95
Fax: +31-(0)-36 521 99 90
E-mail: export@mercado.se