



THE REAL-CHAIRS ARE CE-CERTIFIED
IN ACCORDANCE WITH THE MPA STATUTES
OF MEDICAL DEVICES, LVFS 2003:11.

⊕ REAL 8200 PLUS EL MAMMOGRAPHY

Comfortable, ergonomic for the user and the operator



A mammography exam can be an emotional experience for the patient. That is why it is important to use the right equipment that will make the exam as effective and as pleasant as possible.

Some patients need to stay seated through the entire consultation. It is also easier to get the images right if the patient is seated in a stable, ergonomically correct, seated position. For more time consuming examinations the patient should be seated comfortably and relaxed.

REAL 8200 PLUS EL MAMMOGRAPHY is an examination chair that offers the patient a correct posture and at the same time facilitates a good work environment for the staff. It also makes the examinations quicker and more cost effective. The chair is tried and tested in real environments, developed in collaboration with our global customers and is produced in Sweden.

Base with footbrake and 1/4 foldable foot ring.



Push handle and wired remote control.

⊕ REAL 8200 PLUS EL MAMMOGRAPHY 24V 1G

*Comfortable and ergonomic
for the patient*

- ErgoMedic seating system
- High backrest for extra support
- Padded, adjustable armrests
- Foldable foot ring
- Solid design

*...and facilitates a good work
environment for the staff.*

- Electrical functions for easy adjustment
- 28 cm height adjustment (47-75 cm)
- Large, low friction castors
- Push handle
- Foot brake
- Hygienic upholstery
- Maintenance free
- Makes the examinations more efficient



NEUTRAL POSITION



ELECTRIC SEAT TILT



ELECTRIC BACKREST TILT

REAL 8200 PLUS EL Standard/Screening Mammography



SPECIFICATIONS:

Seat ErgoMedic 44x43 cm Skai Lavender
Without seat tilt
High backrest ErgoMedic 40x34 cm Skai Lavender
Fixed backrest mechanism w. push handle
Armrest 25 cm Skai Black (alt. 30 cm)
High seat height el. liftomat 45-73 cm (alt. low 37-53 cm)
Wired remote control
Base, PLUS clinic white 48x53 cm w. foot brake
Castor 100 mm XL (alt. 100 mm SL low friction)
Brake lever 35 cm
¼ foldable foot ring

Art.nr: 82100901PF

REAL 8200 PLUS EL Advanced/Clinical Mammography



SPECIFICATIONS:

Seat ErgoMedic 44x43 cm Skai Lavender
Electric seat tilt 15° forward/15° backwards
Extra high backrest ErgoMedic 42x43 cm Skai Lavender
Electric backrest mechanism w. push handle
Armrest 25 cm Skai black (alt. 30 cm)
High seat height el. liftomat 47-75 cm (alt. låg 42-62 cm)
Wired remote control
Base, PLUS clinic white 48x53 cm w. foot brake
Castor 100 mm XL (alt. 100 mm SL low friction)
Brake lever 35 cm
Necksupport ErgoMedic PLUS 30x8 cm Skai Lavender
¼ foldable foot ring

Art.nr: 82100902PF

Mercado Medic Sweden

Mercado Medic AB
Tryffelslingan 14 | S-181 57 Lidingö
Postadress: Box 1074, 181 22 Lidingö

Phone: +8-555 143 00
www.mercado.se/en/
kundtjanst@mercado.se

Mercado Medic Nederland

Mercado Medic Nederland
BV Nobelweg 22 3899 BN Zeewolde

Phone: +31 36 521 99 95
www.mercado-medic.nl
admin@mercado-medic.nl

Mercado Medic England

Mercado Medic England
Units 1 & 2, WIndsor Industrial Estate
Hawkins Lane DE14 1QF Burton on Trent

Phone: +44 844 561 6874
www.mercado-medic.co.uk
info@mercado-medic.co.uk

Mercado Medic Germany

Mercado Medic Deutschland GmbH
Dieselstraße 12a , 50259 Pulheim

Phone: +49 223 884 48 640
Mobile: +49 176 806 77 626
www.mercadomedic.de
info@mercadomedic.de



Mercado Medic AB
Tryffelslingan 14
181 57 Lidingö
Sweden

Phone & Fax
P + 46 8 555 143 00

Email & Web
E info@mercado.se
W www.mercado.se