

Category 1 – 6

Below you will find explanation of what the different categories implies.

Category 1

Seat / Leg divider / Side support

- Size: Standard sizes on seats
- Upholstery: Standard shape on foam without changes
- Wooden plate: Standard shape and size
- Fabrics: Within MPAB collection (Outside of standard combination)

Backrest / Trunk support / Neckrest

- Size: Standard size on backrest
- Upholstery: Standard shape on foam without changes
- Wooden plate: Standard shape and size
- Fabrics: Within MPAB collection (Outside of standard combination)

Category 2

Seat / Leg divider / Side support

- Size: Standard sizes on seats up to max 48x53cm
- Upholstery: Build up/hollowed out or in any other way adjusted new shape
- Wooden plate: Cut/drilled or in any other way adjusted shape
- Fabrics: Outside of MPAB collection

Backrest / Trunk support / Neckrest

- Size: Standard size on backrest up to max 46x43cm
- Upholstery: Build up/hollowed out or in any other way adjusted new shape
- Wooden plate: Cut/drilled or in any other way adjusted shape
- Fabrics: Outside of MPAB collection

Category 3

Seat / Leg divider / Side support

- Size: New size up to 56x53cm
- Upholstery: Standard shape on upholstery
- Wooden plate: Standard shape (up to 56x53cm for seats)
- Fabrics: Within MPAB collection (Outside of standard combination)

Backrest / Trunk support / Neckrest

- Size: New size (up to 50x60cm for backrests)
- Upholstery: Standard shape on upholstery
- Wooden plate: Standard shape (up to 50x60cm for backrests)
- Fabrics: Within MPAB collection (Outside of standard combination)

Category 4

Seat / Leg divider / Side support

- Size: New size (up to 56x53cm for seats)
- Upholstery: Build up/hollowed out or in any other way adjusted new shape
- Wooden plate: Cut/drilled or in any other way adjusted shape
- Fabrics: Outside of MPAB collection

Backrest / Trunk support / Neckrest

- Size: New size (up to 50x60cm for Backrest)
- Upholstery: Build up/hollowed out or in any other way adjusted new shape
- Wooden plate: Cut/drilled or in any other way adjusted shape
- Fabrics: Outside of MPAB collection

Category 5

Unique article

- Mechanical components that maintain the products functionality, stability or performance like a custom-made accessory.

Category 6

Custom-made device

- Mechanical components that change the base products function, stability or performance.
- Special combination of standard article that changes the base products function, stability, or performance.
- User characteristics: Movements, non-uniform loads or different seating that demands reinforcements or stabilizing actions needed for the specific product.
- Upholstery: Higher than 12cm. Seats: Size outside of 56x53cm. Backrests: Sizes outside of 50x60cm. These sizes normally need reinforced mechanical parts.