



ECN - Tekniskt ändringsmeddelande

ECN-S08910

Datum: 2024-10-02

Utfärdare: Daniel Salazar

Sida 1 av 2

NYTT UTFÖRANDE AV FÖRGRENINGSKABEL TILL REAL 9200 TWIN EL 24V

Information om ändringar i Mercado Medics produkter eller tjänster.



Bild 1 vänster



Bild 2 höger

Sammanfattning

Förgreningskabeln har fått nytt utförande, men behåller samma artikelnummer.

Artikelnummer som påverkas

805279 Twin Förgreningskabel Kontrollbox-Batteri

Bakgrund

Tillverkningsmetoden har ändrats för att säkerställa tillgängligheten.

Detaljer kring ändringen

Bild 1 till vänster, sida 1 har krympslang runt förgrening enligt nytt utförande. Bild 2 till höger, sida 1 har en formgjuten förgrening som utgår.

Planerad tidpunkt för införande

Q1 2025 med reservation för slutanvändning av artiklar i lager.

Kompatibilitet

Ändringen påverkar inte kompatibiliteten.

Påverkan på kända risker

Ingen påverkan på kända risker.

Frågor

Vid frågor, vänligen kontakta vår produktchef Ronny Fogelqvist.

Kontaktuppgifter:

+46 708 27 96 14

ronny.fogelqvist@mercado.se



ECN - Engineering Change Notice

ECN-S08910

Date: 2024-10-02

Issued by: Daniel Salazar

Page 1 of 2

NEW EXECUTION OF SPLIT CABLE TO REAL 9200 TWIN EL 24V

Information about changes in Mercado Medic's products or services.



Image 1 left



Image 2 right

Summary

The split cable has been redesigned, but retains the same article number.

Affected article numbers

805279

Twin Cable Control Box-Battery

Background

The manufacturing method has been changed to ensure availability.

Description of the change

Image 1 on the left, page 1 has shrink tubing around the branch according to the new design.

Image 2 on the right, side 1 has a molded branch that expires.

Estimated time for implementation

Q1 2025 with reservation for end use of items in stock.

Compatibility

The change does not affect compatibility.

Safety information

No impact on known risks.

Questions

If you have any questions, please contact our Product Manager Ronny Fogelqvist.

Contact Information:

+46 708 27 96 14

ronny.fogelqvist@mercado.se