



Restrictions on CE Marking Regarding Third-Party Leg supports

The limitations specified below apply to the use of third-party leg supports while maintaining the CE marking for the product series REAL 9000 PLUS (excluding REAL 9700/9800 PLUS), REAL 9200 TWIN, REAL 9500/9600 PLUS, REAL 6100 PLUS, and REAL 6100 PLUS LiNX from Mercado Medic AB.

- Third-party leg supports and amputation supports must fit into Mercado Medic's original leg support mount on REAL chairs. Alternatively, they must fit into an adapter manufactured by Mercado Medic. If a specific type of intermediate adapter not currently available is needed, it can be ordered via Mercado Custom Request.
- Manual leg supports may be installed.
- If third-party electric leg supports are required on the electric wheelchair models REAL 6100 PLUS or REAL 6100 PLUS LiNX, beyond those directly marketed by Mercado Medic, can also be ordered via Mercado Medic's Custom Request.

If the above criteria are met, the CE marking can be maintained on the product. This means that the product can be restored to its original condition without requiring Mercado Medic's approval.

However, the product still has a maximum expected service life of 10 years, after which it requires an inspection according to chapter 6.3 in the Instructions for use for the above-mentioned products.