

Restrictions within CE marking regarding Third-Party Seat and Back cushions

The restrictions specified below apply to the use of third-party seat and back cushions while maintaining the CE marking on the following product series from Mercado Medic AB; REAL 9000 PLUS, REAL 9200 TWIN, REAL 6100 PLUS and REAL 6100 PLUS LiNX.

Seat cushions

- External dimensions for REAL 9000 and REAL 6100: maximum width 48 cm, maximum depth 53 cm and maximum height 12 cm.
External dimensions for REAL 9200: maximum width 72 cm and maximum depth 50 cm.
- Upholstered components must meet applicable flammability and environmental requirements.
- Only wooden plates and anti-slip wooden plates supplied by Mercado Medic AB may be used as the base surface.
- It is possible to order wider wooden plates up to 55 cm in width from Mercado Medic as a Unique Article, which still allows CE marking to be maintained.
- Products require Velcro fastening between the wooden plate and the seat cushion. This must be installed according to the cushion supplier's instructions.

Back cushions

- External dimensions of the back cushion may be a maximum width of 50 cm and maximum height of 60 cm.
- Upholstered components must meet applicable flammability and environmental requirements.
- Only wooden back plates and anti-slip back plates (Unique Article) from Mercado Medic AB may be used as the base surface.
- Products require Velcro fastening between the back plate and the back cushion. This must be installed according to the cushion supplier's instructions.
- It is possible to order unique wooden plates up to 50 cm width and 60 cm height from Mercado Medic as a Unique Article.
- Complete third-party back systems, including the supporting structure, may be used following inspection and approval by Mercado Medic AB.
- Adjustable (wheelchair-type) back cushions with clamp brackets, must be mounted on Mercado Medic's Slingback Back Mechanism. There are four models available:
804847 = 39×30 cm, 804849 = 39×40 cm, 804848 = 45×30 cm, 804850 = 45×40 cm
If other sizes are required, they can be ordered as a Unique Article.

If the above criteria are met, the CE marking can be maintained on the product. This means that the product can be restored to its original condition without requiring Mercado Medic's approval.

However, the product still has a maximum expected service life of 10 years, after which it requires an inspection according to chapter 6.3 in the Instructions for use for the above-mentioned products.