



Restrictions within CE marking regarding Third-party belts, harnesses and vests

The restrictions specified below apply to the use of third-party belts, harnesses and vests while maintaining the CE marking on the following product series from Mercado Medic AB; REAL 9000 PLUS, REAL 6100 PLUS and REAL 6100 PLUS LiNX.

- Mounting with Mercado Medic's brackets and screws under the seat for belts with up to 30 mm wide fasteners.
- In the case of wider fastening straps than 30 mm, the bracket must have a minimum of 3 screws in the wooden plate.
- For mounting with Mercado Medic's track system, 1 screw per bracket is sufficient for waist-mounted 2- and 4-point fastening.
- For vests and H-harnesses, upper attachment, Mercado Medic's back track system is required in combination with upper attachment for H-harness (Article number 804184 or complete kit HB2540).

If the above criteria are met, the CE marking can be maintained on the product. This means that the product can be restored to its original condition without requiring Mercado Medic's approval.

However, the product still has a maximum expected service life of 10 years, after which it requires an inspection according to chapter 6.3 in the Instructions for use for the above-mentioned products.