ECN - Tekniskt ändringsmeddelande
ECN-S14593
Datum: 2026-09-02
Utfärdare: Marco Alonso Östberg

REVIDERAT HÖFTBÄLTESFÄSTE

Information om ändringar i Mercado Medics produkter eller tjänster. ECR(s): S14593



Sammanfattning

Fästet för höftbälte (Höftbälte fäste, art. 800878) har reviderats för att även kunna monteras i sidostödsskena.
Bredd och höjd har ökats något samt ett extra hål för montering i skena.

Artikelnummer som påverkas

800878	Höftbälte fäste (artikelnummer oförändrat)

Bakgrund

Behov av montering i sidostödsskena.

Planerad tidpunkt för införande

Under december månad 2026.

Kompatibilitet

Påverkas inte.

Påverkan på kända risker

Inga kända risker påverkas.

Frågor

Vid frågor, vänligen kontakta vår produktchef Ronny Fogelqvist.

Kontaktuppgifter:

+46 708 27 96 14
ronny.fogelqvist@mercado.se

ECN - Engineering Change Notice
ECN-S14593
Date: 2026-09-02
Issued by: Marco Alonso Östberg

REVISED HIP BELT BRACKET

Information about changes in Mercado Medic's products or services. ECR(s):S14593



Summary

The hip belt bracket (item no. 800878) has been revised so that it can also be mounted in the side support rail. The width and height have been increased slightly, and an additional hole has been added for mounting in the rail.

Affected article numbers

800878 Hip belt bracket (article number unchanged)

Background

A need to mount the bracket in the side support rail.

Description of the change

The hip belt bracket has been revised to allow mounting in the side support rail.

Estimated time for implementation

During December 2026.

Compatibility

Not affected.

Safety information

No known risks are affected.

Questions

If you have any questions, please contact our Product Manager Ronny Fogelqvist.

Contact Information:

+46 708 27 96 14
ronny.fogelqvist@mercado.se